



产品名称： MCM-22

制备日期： 2018 年 8 月 3 日

Items	Unit	Data
SiO ₂ /Al ₂ O ₃	mol/mol	30
Na ₂ O Content	wt%	0.02
Crystal Size	μm	0.2-0.5
BET Surface Area	m ² /g	600
Pore Size	nm	0.8
Bulk density	Kg/m ³	500
X-ray Crystallinity	%	90
Pore Capacity	ml/g	>0.30
Lattice constant (nm)	Å	/
Ignition Loss (550°C 3h)	wt%	<8
Shape	D=2-2.5mm, L=5-10mm	
Crush Strengsh	N	>90N

Nanjing JCNano Technology Co., Ltd

Add: Nnijing Road Jianning No. 241 Room 703


Zip: 210015

Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net

XRD 图:



Nanjing JCNano Technology Co., Ltd

Add: Nanjing Road Jianning No. 241 Room 703



Zip: 210015

Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net

SEM:



Nanjing JCNano Technology Co., Ltd

Add: Nnijing Road Jianning No. 241 Room 703

Zip: 210015

Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net



JCNO
吉倉納米

Nanjing JCNano Technology Co., Ltd

Add: Nanjing Road Jianning No. 241 Room 703

Zip: 210015

Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net