




羧基化石墨烯 JCG-1-3-COOH 表征

羧基化石墨烯 JCG-1-3-COOH 参数见下表:

Carboxyl Graphene	
Diameter	1~5 μm
Thickness	0.8~1.2 nm
Single layer ratio	>80%
Carboxyl ratio	5.0 wt%

TEM 投射电镜图



Nanjing JCNano Technology Co., Ltd

Add: Nanjing Road Jianning No. 241 Room 703

Zip: 210015


Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net



红外线图:



Nanjing JCNano Technology Co., Ltd

Add: Nnijing Road Jianning No. 241 Room 703

Zip: 210015

Tel: +86 25 58830538, 400-025-7300

Fax: +86 25 83140772

www.jcno.net